



registered by using the registration code in Reagent Registration Card within the reagent kit.

#### Calibration

Calibration must be performed with the reagents and calibrators from Manufacturer.

Each bottle of lyophilized calibrator should be re-dissolved by distilled water, deionized water, or purified water to 1.0mL before use. Mix it thoroughly and avoid foam formation. Keep dissolved calibrators upright for about 15 minutes to equilibrate to room temperature (15-30°C).

Calibration is necessary under the following conditions:

- (1). New reagent lot number applied,
- (2). **Same lot of reagents used more than 4 weeks,**
- (3). QC out of control,
- (4). When necessary, for example: system maintenance or repair that may affect the analytic performance.

#### Test Procedure

1. Remove the kit from storage at 2-8℃ and take out the required number of reagent cartridges. Allow the reagent cartridges, specimens, calibrators, and/or controls to equilibrate to room temperature (15-30℃) prior to testing.
2. Use each reagent cartridge for each sample, calibrator or control.
3. Place the reagent strips into the reagent rack.
4. Manually pipette the Samples, Calibrators, Controls into the sample well of reagent strip.

**Please note that at least 75µl sample has to be pipetted manually into the reagent cartridge and 50µl sample at least sample is needed for each test.**

***\* Before pipetting: Mix the samples using a vortex-type mixer if necessary, and ensure that the samples and sample diluent (if applicable) DO NOT have any bubbles. The same for calibrators and controls.***

5. Put the sample rack into the reagent Chamber of the ACCRE analyzer.
6. Place the assay tips into the reagent Chamber of the ACCRE analyzer.
7. Input Sample ID, sample type, dilution etc. into the software as directed by USER MANUAL.
8. The Cortisol (CLIA) assay will be performed automatically by the analyzer.
9. After the assay is completed, remove the reagent rack from the instrument. Dispose the used reagent strip and assay tips into an appropriate recipient.

#### Quality control

Each bottle of lyophilized quality control should be re-dissolved by distilled, deionized, or purified water to 1.0mL before use. Mix it thoroughly and avoid foam formation. Keep dissolved quality controls upright for about 15 minutes to equilibrate to room temperature (15-30°C). Quality control can be performed in accordance with local regulations or requirements related to accreditation, as well as requirements defined in the laboratory's quality control procedure.

#### [RESULTS AND INTERPRETATION]

##### Calculation of Results

The results are automatically calculated by the ACCRE analyzer using stored calibration curve which generate by 2-point calibration, and calculated according to a predefined mathematical model, the results expressed in µg/dL.

##### Interpretation of Results

The detection range is from 0.054 µg/dL to 63.400 µg/dL. For samples with

Cortisol below the upper limit of detection, quantitative determination can be performed; if the sample is lower than the lower limit of detection, the reported result is <0.054 µg/dL; if the sample is higher than the upper limit of detection, the reported result is >63.400 µg/dL. Results of Cortisol (CLIA) assays were determined by testing serum specimens from a study in clinical centers with group of totally 566 healthy individuals. The 97.5th percentile reference interval of normal tested value is 4.880 µg/dL to 19.655 µg/dL from 6-10 a.m. and 2.432 µg/dL to 11.934 µg/dL from 4-8 p.m.

***\* Please note that results may differ between laboratories due to variations in population and test method. If necessary, each laboratory should establish its own reference range.***

In interpreting the results, patient's clinical condition including symptoms, diseases history, and other relevant data and information should be considered.

#### [LIMITATIONS]

- 1.For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other tests, clinical impressions, etc.
2. Please note that heterophile antibody in human serum and plasma can interfere the test result.
3. High concentration of HAMA can also cause a false positive or false negative result.
4. If the test results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

#### [PERFORMANCE CHARACTERISTICS]

##### Detection Limit

LoB (Limit of Blank)= 0.036 µg/dL

The Limit of Blank was determined in accordance with the CLSI® EP17-A2 recommendations.

##### Measurement range

The measurement range is the range of values corresponding to the acceptable performance limits (precision and linearity).

**The measurement range of the Cortisol (CLIA) assay is 0.054 µ g/dL to 63.400 µg/dL.**

##### Linearity

Linearity was evaluated according to CLSI ® EP06-A recommendations. The Cortisol (CLIA) assay is between 0.054 µ g/dL ~ 63.400 µ g/dL. The absolute linear correlation coefficient (r) should be bigger than 0.9900.

##### Precision

The precision study was performed according to CLSI® EP05-A3 recommendations. The within-lot precision CV% is less than 10%.

##### Interfering Substances

Potential interference substances were studied according to CLSI ® EP07-A3 recommendations. No significant interference was detected up to the maximum concentrations tested.

| Substance         | Concentration |
|-------------------|---------------|
| Hemoglobin        | ≤ 200 mg/dL   |
| Bilirubin         | ≤ 20 mg/dL    |
| Triglyceride      | ≤ 1000 mg/dL  |
| Rheumatoid factor | ≤ 200 IU/mL   |

#### Accuracy

The relative deviation between the measured value and the calibration concentration is within the range of ±15%.

#### HAMA

Patient samples containing human anti-mouse antibodies (HAMA) may give falsely elevated or decreased values. Although HAMA-neutralizing agents are added, extremely high HAMA serum concentrations may occasionally influence results.

#### [WASTE DISPOSAL]

Dispose used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

#### [INDEX OF SYMBOLS]

| Symbol | Meaning                                                             | Symbol | Meaning                                                           | Symbol | Meaning                                       |
|--------|---------------------------------------------------------------------|--------|-------------------------------------------------------------------|--------|-----------------------------------------------|
|        | Consult instructions for use                                        |        | Contains sufficient for <n> tests                                 |        | Date of manufacture                           |
|        | In vitro diagnostic medical device                                  |        | Use-by date                                                       |        | Do not re-use                                 |
|        | Temperature limit                                                   |        | Batch code                                                        |        | Catalogue number                              |
|        | Keep away from sunlight                                             |        | Keep dry                                                          |        | Manufacturer                                  |
|        | Authorized representative in the European Community/ European Union |        | Do not use if package is damaged and consult instructions for use |        | Contains biological material of animal origin |
|        | Biological risks                                                    |        | Caution                                                           |        | CE marking                                    |
|        | Upright                                                             |        | Fragile, handle with care                                         |        | Do not roll                                   |
|        | Control                                                             |        |                                                                   |        |                                               |

#### [REFERENCES]

1. Turpeinen U, Hämäläinen E. Determination of cortisol in serum, saliva and urine[J]. Best Practice & Research Clinical Endocrinology & Metabolism, 2013, 27(6): 795-801.
2. Beisel W R, DiRaimondo V C, Forsham P H. Cortisol Transport and Disappearance[J]. Annals of Internal Medicine, 1963, 58(4): 722-723.
3. Gatti R, Antonelli G, Prearo M, et al. Cortisol assays and diagnostic laboratory procedures in human biological fluids[J]. Clinical biochemistry, 2009, 42(12): 1205-1217.
4. Moore A, Aitken R, Burke C, et al. Cortisol assays: guidelines for the provision of a clinical biochemistry service[J]. Annals of clinical biochemistry, 1985, 22(5): 435-454.
5. Tsigos C, Chrousos G P. Hypothalamic-pituitary-adrenal axis, neuroendocrine factors and stress[J]. Journal of psychosomatic research, 2002, 53(4): 865-871.

6. Boscatto L M, Stuart M C. Heterophilic antibodies: a problem for all immunoassays[J]. Clinical chemistry, 1988, 34(1): 27-33.
7. Bjerner J, Nustad K, Norum L F, et al. Immunometric assay interference: incidence and prevention[J]. Clinical chemistry, 2002, 48(4): 613-621.
- 8.Boscatto Lm,Stuart MC.Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.
- 9.Kricka L.Interferences in immunoassays - still a threat. Clin Chem 2000; 46:1037-1038.
- 10.Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002;48:613-621.



Shenzhen Tisenc Medical Devices Co., Ltd.  
Address: 11F, 11G, Kechuang Building, Quanzhi Science and Innovation Park, Maozhoushan Industrial Park, Shajing Community, BaoAn District, Shenzhen, 518104, P.R.China  
Phone: +86-400-830-8768 +86-0755-23225549  
E-mail: sales@tisenc.com



CMC Medical Devices&Drugs S.L.  
C/HoracioLengono18,CP29006,Málaga,Spain

Eng 20220402 -Ver 2.0