



# ichroma™ Anti-HBs

## INTENDED USE

**ichroma™ Anti-HBs** is a fluorescence Immunoassay (FIA) for the qualitative determination of antibody to Hepatitis B surface antigen (anti-HBs) in human whole blood/serum/plasma. It is useful as an aid to diagnosis of HBV infection or following Hepatitis B virus vaccination.

For *in vitro* diagnostic use only.

## INTRODUCTION

Viral hepatitis is a serious global health problem affecting over two billion people worldwide and approximately one million people die each year due to cirrhosis of the liver and hepatocellular carcinoma (HCC), which are commonly associated with chronic hepatitis. The majority of hepatitis viral infections are caused by three distinct virus types: Hepatitis A (HAV, Hepatitis B(HBV), Hepatitis C (HCV).<sup>1,2</sup> The risk of developing chronic infection by HBV varies inversely with age and is highest for infants infected at birth compared to older children and adults. Up to 90% of infants infected with HBV will develop chronic infection leading to cirrhosis of the liver or HCC compared to 6-10% of adults who acquire HBV infection.<sup>3</sup> Determination of antibodies directed against HBV surface antigen (anti-HBs) is used to evaluate a person's immune status to HBV infection or to aid in the laboratory diagnosis of HBV infection when used in conjunction with other laboratory methods. The test is performed to assess the need for vaccination (if anti-HBs is absent or below levels considered protective), following completion of vaccination against HBV in high risk groups (healthcare workers, chronic renal failure patients, HIV infected persons), or to monitor recovery from acute HBV infection. The presence of anti-HBs following acute infection generally indicates recovery and immunity from infection.

## PRINCIPLE

The test uses a sandwich immunodetection method.

The detector antigen in buffer binds to antibody in sample, forming antigen-antibody complexes, and migrates onto nitrocellulose matrix to be captured by the other immobilized antigen on test strip.

More antibodies in the sample will form the more antigen-antibody complexes which lead to stronger fluorescence signal by detector antigen, which is processed by the instrument for ichroma™ tests to display the 'Positive' / 'Negative' / 'Indeterminate' in the sample.

## COMPONENTS

**ichroma™ Anti-HBs** consists of 'cartridges', 'detector tubes' and 'detector diluent'.

- The cartridge contains the membrane called a test strip which has recombinant HBsAg at the test line and chicken IgY at the control line. All cartridges are individually sealed in an aluminum foil pouch containing a desiccant in a box.
- The detector tube has a granule containing recombinant HBsAg-fluorescence conjugate, anti-chicken IgY-fluorescence conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide as a preservative in Tris-HCl. All detector tubes are packed in a box.
- The detector diluent contains salt and bovine serum albumin (BSA) as a stabilizer and sodium azide as a preservative in Tris-HCl and it is pre-dispensed in a vial. The detector diluent is packed in a box.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instructions for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge, detector tube, detector diluent and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges or detector tubes. A cartridge should be used for testing one sample only. A detector tube should be used for processing of one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use the cartridge, if pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- If test components and/or sample are stored in refrigerator, then allow cartridge, detector tube, detector diluent and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for ichroma™ tests may generate slight vibration during use.
- Used cartridges, detector tubes, detector diluent and pipette tips should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- The detector tube and the detector diluent contain sodium azide (NaN<sub>3</sub>), and it may cause certain health issues like convulsions, low blood pressure, low heart rate, loss of consciousness, lung injury and respiratory failure.
- **ichroma™ Anti-HBs** will provide accurate and reliable results subject to the below conditions.
  - **ichroma™ Anti-HBs** should be used only in conjunction with the instrument for ichroma™ tests.
  - Have to use recommended anticoagulant.

Recommended anticoagulant

Na EDTA, K<sub>2</sub> EDTA, K<sub>3</sub> EDTA,  
Na-heparin, Li-heparin, Sodium citrate

## STORAGE AND STABILITY

| Storage condition |                     |            |            |
|-------------------|---------------------|------------|------------|
| Component         | Storage Temperature | Shelf life | Note       |
| Cartridge         | 4 – 30 °C           | 20 months  | Disposable |
| Detector tube     | 2 – 8 °C            | 20 months  | Disposable |
| Detector diluent  | 2 – 8 °C            | 20 months  | Unopened   |
|                   |                     | 20 months  | Opened     |

- After the cartridge pouch is opened, the test should be performed immediately.

## LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigens to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigens with time and/or temperature may also cause false negative result as it makes the antigens unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

## MATERIALS SUPPLIED

### REF CFPC-52

Components of **ichroma™ Anti-HBs**

- Cartridge Box:
  - Cartridge 25
  - Detector tube (packaged in a separated box) 25
  - Detector diluent (packaged in a separated box) 1
  - ID chip 1
  - Instructions for use 1

## MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **ichroma™ Anti-HBs**.

Please contact our sales division for more information.

- Instrument for **ichroma™** tests
  - ichroma™ II** REF FPRR021
  - ichroma™ III** REF FPRR037
  - ichroma™ M2** REF FPRR031
- Boditech Anti-HBs Control** REF CFPO-144

## SAMPLE COLLECTION AND PROCESSING

The sample type for **ichroma™ Anti-HBs** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection. when collected sample is stored at room temperature.

- The samples (serum, plasma) should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.
- The samples (whole blood, serum, plasma) may be stored for a week at 2-8 °C prior to being tested. If testing will be delayed more than a week, samples (serum, plasma) should be frozen at -20 °C.
- The samples (serum, plasma) stored frozen at -20 °C for 2 months showed no performance difference.
- However, the whole blood sample should not be kept in a freezer in any case.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.

## TEST SETUP

- ichroma™ Anti-HBs** should be performed at environment temperature 25 ± 3 °C.
  - Check the contents of **ichroma™ Anti-HBs**: Sealed cartridges, detector tubes, detect diluent, ID chip and Instructions for use.
  - Ensure that the lot number of the cartridge matches that of the detector tube, the detector diluent as well as an the ID chip.
  - If the sealed cartridge, the detector tube and the detector diluent have been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
  - Turn on the Instrument for **ichroma™** tests.
  - Insert the ID chip into the 'ID chip port'.
- ※ Please refer to the instrument for **ichroma™** tests operation manual for complete information and operating instructions.

## TEST PROCEDURE

### ► **ichroma™ II, ichroma™ M2**

#### **Multi test mode / Read now mode**

- Select the sample type (whole blood or serum/plasma) on the screen.
- Take 100 µL of detector diluent using a pipette and dispense it to the detector tube containing a granule. When the granule form is completely dissolved in the tube, it becomes detection buffer.
- Take 50 µL of sample (whole blood/serum/plasma/control) using a pipette and dispense it to the detector tube.
- Close the lid of the detector tube and mix the sample thoroughly by shaking it about 10 - 20 times. (The sample mixture must be used immediately.)
- Take 100 µL of the sample mixture and dispense it into the sample well of the cartridge.
- Leave the cartridge at room temperature for 15 minutes.

△ Scan the sample loaded cartridge immediately when the incubation time is over. If not, it will cause inaccurate test result.

- To scan the sample-loaded cartridge, insert it into the cartridge holder of the instrument for **ichroma™** tests. Ensure proper orientation of the cartridge before pushing it all the way inside the cartridge holder. An arrow is marked on the cartridge especially for this purpose.

- 8) Press the 'Select' or tap the 'Start' button on the instrument for ichroma™ tests to start the scanning process.  
(ichroma™ M2 will start the test automatically after inserting.)
- 9) The instrument for ichroma™ tests will start scanning the sample loaded cartridge immediately.
- 10) Read the test result on the display screen of the instrument for ichroma™ tests.

#### Single test mode/ Walk away mode

- 1) The test procedure is same with the 'Multi test mode 1) – 5)'
- 2) Insert the sample-loaded cartridge into the holder of the instrument for ichroma™ tests. Ensure proper orientation of the cartridge before pushing it all the way inside the cartridge holder. An arrow is marked on the cartridge especially for this purpose.
- 3) Press the 'Select' or tap the 'Start' button on the instrument for ichroma™ tests.  
(ichroma™ M2 will start the test automatically after inserting.)
- 4) The cartridge goes inside the instrument for ichroma™ tests and will automatically start scanning the sample-loaded cartridge after 15 minutes.
- 5) Read the test result on the display screen of the instrument for ichroma™ tests.

#### ► ichroma™ III

- 1) The test procedure is same with the 'Single test mode'.

#### - Cross-reactivity

There was no significant cross-reactions in the **ichroma™ Anti-HBs** test measurement.

| No. | Cross-reactivity materials  |
|-----|-----------------------------|
| 1   | Cytomegalovirus (CMV)       |
| 2   | Epstein-Barr virus (EBV)    |
| 3   | Hepatitis A virus(HAV)      |
| 4   | Hepatitis C virus(HCV)      |
| 5   | Herpes simplex virus(HSV)   |
| 6   | Rubella virus               |
| 7   | Varicella-zoster virus(VZV) |
| 8   | Treponema pallidum          |
| 9   | Anti Nuclear Antibody(ANA)  |
| 10  | Rheumatoid factor(RF)       |
| 11  | Early stage of pregnancy    |
| 12  | Middle stage of pregnancy   |

#### - Interference

There was no significant interference from these material with the **ichroma™ Anti-HBs** test.

| No. | Interference materials | Conc.       |
|-----|------------------------|-------------|
| 1   | Heparin                | 100,000 U/L |
| 2   | EDTA                   | 4 µM        |
| 3   | Sodium citrate         | 0.085 M     |
| 4   | Hemoglobin             | 2 mg/mL     |
| 5   | BSA                    | 60 mg/mL    |
| 6   | Bilirubin              | 400 µM      |
| 7   | Triglycerides          | 1.5 mg/mL   |
| 8   | Cholesterol            | 20 mM       |

#### INTERPRETATION OF TEST RESULT

- The instrument for ichroma™ tests calculates the test result automatically and displays "Positive" / "Negative" / "Indeterminate".

| Result (mIU/mL) | Result                    | Note                           |
|-----------------|---------------------------|--------------------------------|
| ≤ 5             | Negative<br>for anti-HBs. | No need<br>to additional test. |
| 5 < Titer < 15  | Indeterminate.            | Need to retest.                |
| ≥ 15            | Positive<br>for anti-HBs. | Need<br>to confirmation test.  |

#### QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **ichroma™ Anti-HBs**. For more information regarding obtaining the control materials, contact [Boditech Med Inc.'s Sales Division for assistance](#).  
(Please refer to the instructions for use of control material.)

#### PERFORMANCE CHARACTERISTICS

- **Analytical Sensitivity**
  - Limit of Blank (LOB) 3.60 mIU/mL
  - Limit of Detection (LOD) 4.55 mIU/mL
  - Cut-off 15 mIU/mL
- **Analytical Specificity**

#### ■ Precision

##### - Between lot

One person tested three different lots of **ichroma™ Anti-HBs**, ten times at each concentration of the control standard.

##### - Between person

Three different persons tested same lot of **ichroma™ Anti-HBs**, five times at each concentration of the control standard.

##### - Between day

One person tested one lot of **ichroma™ Anti-HBs** during three days, five times at each concentration of the control standard.

##### - Between site

One person tested one lot of **ichroma™ Anti-HBs** in three different site, five times at each concentration of the control standard.

| Cal      | Between lot                     |                  | Between person                  |                  |
|----------|---------------------------------|------------------|---------------------------------|------------------|
|          | Positive/<br>Number of<br>tests | Positive<br>rate | Positive/<br>Number of<br>tests | Positive<br>rate |
| Negative | 0/30                            | 0%               | 0/15                            | 0%               |
| High     | 30/30                           | 100%             | 15/15                           | 100%             |
| Mid      | 30/30                           | 100%             | 15/15                           | 100%             |
| Low      | 30/30                           | 100%             | 15/15                           | 100%             |
| Cal      | Between day                     |                  | Between site                    |                  |
|          | Positive/<br>Number of<br>tests | Positive<br>rate | Positive/<br>Number of<br>tests | Positive<br>rate |
| Negative | 0/15                            | 0%               | 0/15                            | 0%               |
| High     | 15/15                           | 100%             | 15/15                           | 100%             |
| Mid      | 15/15                           | 100%             | 15/15                           | 100%             |
| Low      | 15/15                           | 100%             | 15/15                           | 100%             |

#### ■ Comparability with reference product

|                      |          | Reference product |          |       |
|----------------------|----------|-------------------|----------|-------|
|                      |          | Positive          | Negative | Total |
| ichroma™<br>Anti-HBs | Positive | 112               | 5        | 117   |
|                      | Negative | 10                | 679      | 689   |
|                      | Total    | 122               | 684      | 806   |

- Positive Comparability: 91.8 %

- Negative Comparability: 99.3 %

- Total Comparability: 98.1 %

#### REFERENCES

1. The Global Burden of Liver Disease: The Major Impact of China. Hepatology. 2014, 60:2099-2108
2. Viral hepatitis in resource-limited countries and access to antiviral therapies: current and future challenges. Future Virol. 2013, 8:371-380
3. Mahoney FJ, et al. Update on diagnosis, management and prevention of hepatitis B virus infection, 1999, Clin Microbiol Rev, 12:351-366.

**Note:** Please refer to the table below to identify various symbols.

|  |                                    |
|--|------------------------------------|
|  | Sufficient for <n> tests           |
|  | Read instruction for use           |
|  | Use by Date                        |
|  | Batch code                         |
|  | Catalog number                     |
|  | Caution                            |
|  | Manufacturer                       |
|  | In vitro diagnostic medical device |
|  | Temperature limit                  |
|  | Do not reuse                       |

For technical assistance, please contact:

**Boditech Med Inc.'s Technical Sales**

Tel: +(82) -33-243-1400

E-mail: TS@boditech.co.kr



**Boditech Med Inc.**

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

Tel: +(82) -33-243-1400

Fax: +(82) -33-243-9373

www.boditech.co.kr

